

Overzicht per hoortoestel Tubes & Receivers Tips & Sleeves



RIC BC LI T 81

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



RIC LI 81

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



RIC LI T 81

- MiniReceiver 3.0 (S,M,P)
- MiniReceiver 3.0 HP
- EarTips 3.0 (Open, Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



IF CIC LI 81

- Sleeves 3.0 (Vented, Closed, Power)



OTE LI X80

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open,Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



OTE X80

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open,Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



P X80

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open,Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



P LI X19

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open,Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



OTE LI X19

- ThinTube 3.0 (Standaard 0,9 mm)
- ThinTube 3.0 (Power 1,4 mm)
- EarTips 3.0 (Open,Tulip)
- Sleeves 3.0 (Vented, Closed, Power)



Overzicht per hoortoestel Tubes & Receivers Tips & Sleeves



RIC LI X19

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Power)
- ClickSleeves (Vented, Closed)



RIC 312 X19

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Power)
- ClickSleeves (Vented, Closed)



IF CIC X19

- ClickSleeves (Vented, Closed)



OTE LI X18

- ThinTube met koppelstukje
- Tips (Open, Closed. Tulip, Double)



OTE X18

- ThinTube met koppelstukje
- Tips (Open, Tulip, Closed, Double)



RIC 312 X18

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Double)
- ClickSleeves (Vented, Closed)



RIC+ 13 X18

- MiniReceiver 2.0 (S,M,P)
- MiniReceiver 2.0 HP
- ClickDomes (Open, Tulip, Double)
- ClickSleeves (Vented, Closed)



MINI RIC 10 X18

- MiniReceiver 2.0 (S,M,P)
- ClickDomes (Open, Tulip, Double)
- ClickSleeves (Vented, Closed)



Open X17

- ThinTube
- EarTips (Open, Closed. Tulip, Double)



Overzicht filters bij de juiste Tubes & Receivers IHO & Molds

WaxGuard 3.0 nanocare

- ThinTube 3.0 (Standaard)
- ThinTube 3.0 (Power)
- MiniReceiver 3.0 (S,M,P)
- Extended earmold 3.0
- IF CIC Li X81

Voor instructie scan of klik op de code



MiniReceiver 3.0

WaxGuard MiniReceiver 2.0

- MiniReceiver 2.0 (S,M,P)
- IF CIC X19



MiniReceiver 2.0

QuickGuard

- ClickMold
- MiniReceiver 2.0 HP
- CIC / IIC X81
- ITE X18



ClickMold

WaxGuard 3.0 XL nanocare cerustop

- MiniReceiver 3.0 HP



MiniReceiver 3.0 HP

MF1 Microfoonfilter

- CIC / IIC model
- in Mokka, Bruin, Donkerbruin en Zwart
- MF1 Removal Tool ook los te bestellen



CIC / IIC X81



Overzicht Tools bij de juiste Tubes & Receivers IHO & Molds


Meettool Earwear 3.0

- Lengte Tube en MiniReceiver 3.0 meten



Connex EarMold 3.0 Measuring card

- Afdruk meten of deze geschikt is voor Maatwerk t.b.v. Earwear 3.0



Earmold 3.0 Removal Tool

- Verwijderen Earmold 3.0 (kort) van een MiniReceiver 3.0
Voor instructie scan of klik op de code



Earmold 3.0 Handling Tool

- Plaatsen en verwijderen van Earmold 3.0 (kort en lang) van een MiniReceiver 3.0
Voor instructie scan of klik op de code



ClickMold Removal Tool

- Verwijderen Clickmold van een MiniReceiver 2.0



MF1 Removal Tool

- Verwijderen en plaatsen MF1 Filter



Tool Removal H2770 BLU

- Verwijderen van de toonbocht bij X19 en X80 hoortoestellen



Overzicht Conchalock 3.0 bij de juiste Tubes & Receivers

ConchaLock 3.0 universeel

- ThinTube 3.0 (Standaard 0,9)
- ThinTube 3.0 (Power 1,4)
- MiniReceiver 3.0 (S,M,P)

Voor instructie scan of klik op de code



Conchalock 3.0 maat S

- ThinTube 3.0 (Standaard 0,9)
- ThinTube 3.0 (Power 1,4)
- MiniReceiver 3.0 S

Voor instructie scan of klik op de code



Conchalock 3.0 maat M

- MiniReceiver 3.0 M



Conchalock 3.0 maat P

- MiniReceiver 3.0 P



Overzicht Tips & Sleeves bij de juiste Tubes en Receivers

EarTips 3.0 (Open, Tulip)

- ThinTube 3.0 (Standaard 0,9)
- ThinTube 3.0 (Power 1,4)
- MiniReceiver 3.0 (S,M,P)

Voor instructie scan of klik op de code



Sleeves 3.0 XS,S,M,L (Vented, Closed, Power)

- ThinTube 3.0 (Standaard 0,9)
- ThinTube 3.0 (Power 1,4)
- MiniReceiver 3.0 (S,M,P)
- IF CIC Li X81

ClickDomes 2.0 (Open, Tulip, Closed, Double)


- MiniReceiver 2.0 (S,M,P)

ClickSleeves 2.0 XS,S,M,L (Vented, Closed)

- MiniReceiver 2.0 (S,M,P)
- IF CIC X19

Tips (Open, Tulip, Closed, Double)

- Open / ThinTube (Open X17)
- Thin Tube met koppelstukje (OTE X18)




Overzicht ThinTubes en miniReceivers

MiniReceivers 3.0

- S Lengte 0 t/m 5
- M Lengte 0 t/m 5
- P Lengte 0 t/m 5
- HP Lengte 0 t/m 5 (via maatwerk bestellen)

Voor instructie scan of klik op de code



ThinTubes 3.0

- Standaard 0,9 mm - Lengte 0 t/m 4
- Power 1,4 mm - Lengte 0 t/m 4



ThinTube 3.0
Standaard 0,9 mm

ThinTube 3.0
Power 1,4 mm

MiniReceivers 2.0

- S Lengte 0 t/m 3
- M Lengte 0 t/m 4
- P Lengte 0 t/m 3
- HP Lengte 0 t/m 3 (via maatwerk bestellen)



ThinTubes

- ThinTube met koppelstukje te bevestigen - Lengte 1 t/m 4
- Open / ThinTube - Lengte 1 t/m 6



ThinTube met koppelstukje

Open X17 / ThinTube

